




# NanoScan OP800

## Nanopositioning Piezo Objective Scanner

The NanoScan OP800 provides fast step and settle times over an 800µm range. Its market leading positioning accuracy and resolution originates from its unique mechanical design and integral capacitive feedback sensors. Compatible with most microscopes and objective lenses, the OP800 has user configurable settings optimized for different objective sizes, weights, and performance needs. You simply select the best setting for your particular application using the software provided.



### Key features

- 800µm closed loop travel range (950µm open loop range typical)
- Rapid settling times even with large mass objectives
- Stiff construction supports heavy objectives with standard settings
- Friction free flexures deliver high stiffness and minimize off-axis motions for the highest repeatability and fastest cycle times
- Tested to perform for more than 10 million full-range cycles
- Capacitive positioning sensors used to give low nanometer positioning and the best resolution, stability (low drift) and repeatability
- Can be used on both upright and inverted microscope
- Connectors with built in stage calibration provide plug and play electronics which can be interchanged, minimizing system down times
- Autocalibration centers the travel range ensuring full range scanning in all appropriate environments

### Applications

- Optical sectioning producing 3D images
- Autofocus systems for time lapse imaging
- High content screening
- Surface analysis
- Wafer inspection
- Scanning interferometry
- Multiphoton Microscopy

## Specifications for NanoScan OP800 Positioner:

Parameter	Specification
<b>Material</b>	Stainless Steel and Aluminum
<b>Closed Loop Range</b>	800 µm
<b>Resolution*</b>	2.5 nm under 250 g load
<b>Linearity*</b>	0.01%
<b>Repeatability*</b>	5 nm Slow PID, 100 µm step, 1SD
<b>Max. Objective Load</b>	500 g (higher loads on request)
<b>Loaded Resonant Frequency 150 g/ 250 g/ 500 g load*</b>	86/82/42 Hz
<b>5% Settle 0.5 µm Step 150 g/250 g/500 g Load*</b>	19 ms/17 ms/21 ms
<b>0.5% (500nm) Settle 100µm Step 150 g/250 g/500 g Load*</b>	51 ms/ 59 ms/70 ms
<b>Cable Lengh</b>	2 m
<b>Objective Clearance</b>	40 mm
<b>Optical Path Length Extension</b>	13 mm; 12 mm with M32

\*Typical specifications

## NanoScan OP800 System Ordering Information

Part Number	Description
<b>QGOP-800-UP-D1</b>	System incorporating NanoScan OP800 (0-500g load) for upright microscopes using objective lenses and NPC-D-6110 controller
<b>QGOP-800-INV-D1</b>	System incorporating NanoScan OP800 (0-500g load) for inverted microscopes using objective lenses and NPC-D-6110 controller

## Microscope Thread Adapter Ordering Information

Part Number	Description
<b>QG-OP-MIC-RMS</b>	Microscope adapter RMS
<b>QG-OP-MIC-M25</b>	Microscope adapter M25 x 0.75
<b>QG-OP-MIC-M32</b>	Microscope adapter M32 x 0.75
<b>QG-OP-MIC-M27</b>	Microscope adapter M27 x 0.75
<b>QG-OP-MIC-W26</b>	Microscope adapter W26 (M26 x 1/36")

## Specifications for NPC-D-6110 Controller

Parameter	Specification
<b>Mechanical</b>	
<b>NPC-D-6110 Controller Box</b>	318 mm x 240 mm 90 mm
<b>Weight of Controller Box</b>	3.0 kg
<b>Cooling</b>	Convection cooled - temperature controlled fans
<b>Electrical</b>	
<b>Power Input</b>	100 to 240 VAC nominal 47 to 63 Hz
<b>Connectivity</b>	
<b>USB - Type B Connector</b>	2.0 compliant
<b>Analog Input Command</b>	BNC - 0 -10V (+/- 10V max)
<b>Analog Position Monitor Output</b>	BNC - 0-10V (+/- 10V max)
<b>"TRIG" input/"TRIG" output "IN-POS" Output and Quadrature Interface</b>	25 pin D-type socket - 5V TTL
<b>Controller Synchronizing Signals</b>	9 pin D-type socket
<b>Enviromental - Operational</b>	
<b>Temperature</b>	10 to 40C
<b>Relative Humidity</b>	5 to 80%

\*Typical specifications

## Objective Lens Adapter Ordering Information

Part Number	Description
<b>QG-OP-OBJ-RMS</b>	Objective adapter M32 x 0.75 to RMS
<b>QG-OP-OBJ-M25</b>	Objective adapter M32 x 0.75 to M25 x 0.75
<b>QG-OP-OBJ-M27</b>	Objective adapter M32 x 0.75 to M27 x 0.75
<b>QG-OP-OBJ-W26</b>	Objective adapter M32 x 0.75 to W26 (M26 x 1/36")

## Parfocal Objective Spacer Ordering Information

Part Number	Description
<b>QG-OP-SPACE-RMS</b>	Objective lens spacer RMS
<b>QG-OP-SPACE-M25</b>	Objective lens spacer M25 x 0.75
<b>QG-OP-SPACE-M32</b>	Objective lens spacer M32 x 0.75
<b>QG-OP-SPACE-M27</b>	Objective lens spacer M27 x 0.75
<b>QG-OP-SPACE-W26</b>	Objective lens spacer W26 (M26 x 1/36")

## NPC-D-6110 Controller




- The powerful digital controller drives the OP800 at the fastest speeds possible.
- Motion control algorithms with acceleration/deceleration control and active damping reduce overshoot.
- Velocity control algorithm gives ultra-smooth ramps for applications such as focus stacking or focus bracketing.
- Market-leading 20  $\mu$ s update rate
- Fastest recovery time between z-stacks, providing enhanced time resolution
- Selectable tuning presets which optimize for step settle with different loads such as incubator


## Interfacing

- Analog command input and position output (0-10 V) for compatibility with existing systems.
- Digital commands over USB for maximum accuracy with a DLL interface for customer software. In-position digital outputs can be used to control camera imaging by providing rapid z-stacking.
- Digital quadrature/step-and-direction commands allow high-speed control with a standard 2-wire motion controller interface, without the need for expensive high-precision ADCs/DACs.
- Playback of custom-programmed waveforms such as constant-velocity profiles. Separate digital trigger outputs can be activated at custom-defined points to control external equipment such as camera imaging.
- Compatible with Queensgate Nanobench, Micro-Manager, Metamorph, Elements via RS232C, and other customer software using the DLL interface provided.
- Can be connected to Prior ProScan III for integrated fine z-control.

## Step and Settle Time



## Dimensional Drawings



电话: 0755-84870203  
邮箱: sales@highlightoptics.com

