



电话: 0755-84870203

邮箱: sales@highlightoptics.com








## Piezo Z-Stage for Märzhäuser Scanning Stages in the SCAN IM Series

Specifications	USB / ANALOG IN¹
Travel range	300 μm
Resolution	4 nm / 18 nm
Repeatability	4 nm / 18 nm (unidirectional)
Linearity	0,5 % / 2,5 %
Settling time (unloaded/under load)	25 ms
Load <sup>2</sup>	max. 500 g
Operation mode	closed-loop operation
Material	aluminium
Weight	1 kg (stage incl. 2 m connecting cable)



## **Order Information**

Order No.: 00-55-550-0800

## **Product Features**

- compatible with nearly all Märzhäuser scanning stages in the SCAN IM series
- direct connection of the piezo amplifier Nano-Drive® C to TANGO controllers with AUX I/O interface
- comfortable operation of your complete positioning system with an only operating device (e. g. joystick)
- high-speed piezo Z-motion with travel ranges up to 300 μm
- direct insertion of well plates without an additional stage insert
- extra large handle recesses for easy handling of well plates with lid
- measuring instrument for easy offset calibration included in delivery

## **Accessories**

Universal stage insert Order No.: 00-55-500-1020

Stage insert for slides

Order No.: 00-55-500-1000 (1 slide) Order No.: 00-55-500-1030 (2 slides)

Stage insert for 1 Petri dish Order No.: 00-55-500-1010



<sup>1</sup> operating mode "USB": connected to PC operating mode "ANALOG IN": connected to TANGO controller <sup>2</sup> recommended value, higher loads on request